

Impact of Percutaneous Closure Device Type On Vascular and Bleeding Complications After TAVR: A Post Hoc Analysis From the BRAVO-3 Randomized Trial

Power, David et al. Catheter and Cardiovasc Interv. 2019 Jun 1;93(7):1374-1381.

PURPOSE

- To investigate the impact of vascular closure devices (VCD) on vascular and bleeding complications after transfemoral TAVR in the large Bivalirudin Versus Heparin Anticoagulation in Transcatheter Aortic Valve Replacement (BRAVO-3) randomized trial population.

METHODS

- A post hoc analysis of the BRAVO-3 trial included stratification of patients based on the type of VCD used. Out of 802 patients, 746 received a VCD:
 - Prostar™ XL: N = 352 (47%)
 - ProGlide™: N = 394 (53%)
- Examined the 30-day incidence of the following:
 - Major or minor vascular complications
 - Major bleeding (BARC \geq 3b)
 - Acute Kidney Injury (AKI)
 - Major Adverse Cardiac and Cerebrovascular Events (MACCE; death, myocardial infarction or stroke)

RESULTS

- The Prostar system was associated with an increased incidence of major or minor vascular complications when compared to the ProGlide system (adjusted OR: 0.54, 95% CI, $p < .01$).
- The overall incidence of major or minor vascular complications was 24% with Prostar and 15% with ProGlide. The incidence of major bleeding (BARC \geq 3b) was 11% and 8% respectively.
- Incidence of Acute Kidney Injury was lower with the ProGlide system as well.
- There was no significant difference between the two VCDs in bleeding, MACCE, and death.

AUTHOR CONCLUSIONS

- Compared to Prostar, the ProGlide VCD was associated with a lower rate of major or minor vascular complications and lower rates of AKI after transfemoral TAVR.
- Operators should be aware of the impact percutaneous VCDs have on vascular and bleeding complications.

Indications for Use. The Early Bird is indicated for the introduction of catheters, catheter balloons, and other diagnostic and interventional devices into the femoral artery or femoral vein while maintaining hemostasis during diagnostic and interventional endovascular procedures.

Contraindications. There are no known contraindications for Early Bird.

Impact of Percutaneous Closure Device Type On Vascular and Bleeding Complications After TAVR: A Post Hoc Analysis From the BRAVO-3 Randomized Trial

(Continued)

30-DAY CLINICAL OUTCOMES

	PROSTAR™ XL N = 352	PROGLIDE™ N = 394
MAJOR OR MINOR VASCULAR COMPLICATIONS	86 (24%)	61 (15%)
MAJOR BLEEDING (BARC ≥ 3B)	39 (11%)	33 (8%)
ANY BARC BLEEDING	167 (47%)	194 (49%)
ACUTE KIDNEY INJURY	70 (25%)	57 (17%)
MACCE	28 (8%)	32 (8%)
DEATH	21 (6%)	16 (4%)

KEY POINTS

- Although there was no significant difference between ProGlide and Prostar, bleeding complications after transfemoral TAVR occurred frequently in patients with either VCD (9.5% major bleeding BARC ≥ 3b).
- Preventative and bleeding avoidance strategies, such as the early detection of bleeding events, can reduce the risk of bleeding in a high-risk TAVR patient population.

Indications for Use. The Early Bird is indicated for the introduction of catheters, catheter balloons, and other diagnostic and interventional devices into the femoral artery or femoral vein while maintaining hemostasis during diagnostic and interventional endovascular procedures.

Contraindications. There are no known contraindications for Early Bird.